

PUTTING ON SILVER WHINNYS®



-1-



-2-

1. If your horse will allow, put a gallon plastic bag over the hoof.
2. Slide the sock over the bag, completely over the hoof and onto the leg.
3. Remove the bag and put the hoof down on the ground.
4. Adjust the sock up and down on the leg. The colored yarn indicates the bottom of the sock. Do not twist the socks. Keep the ribbing straight. PRACTICE MAKES PERFECT!

Use duct tape to cover the rough spots on the hoof if your horse won't tolerate a plastic bag.

ONLY USE TWO SOCKS FOR THE FIRST 24-48 HOURS. GIVE YOUR HORSE TIME TO GET USED TO THEM! THE FIRST TIME YOU PUT THEM ON CAN BE A CHALLENGE UNTIL YOU GET THE HANG OF IT. PRACTICE MAKES PERFECT. YOU SHOULD BE ABLE TO FIT YOUR HAND BETWEEN THE SOCK AND THE CANNON BONE OF THE LEG. THEY ARE NOT MEANT TO BE A COMPRESSION SOCK. IF YOU BELIEVE THEY FEEL TOO TIGHT ON THE LEG, TAKE THE SOCKS OFF AND CALL OR EMAIL US. THIS IS IMPORTANT IF YOUR HORSE HAS SUFFERED CELLULITIS OR LYMPHANGITIS PRESENTLY OR IN THE PAST. ONLY ADD THE SECOND PAIR IF THE FIRST PAIR ARE PERFORMING CORRECTLY.

WE HAVE GREAT TECH SUPPORT. CALL US AT

850-907-5724 OR EMAIL US AT: SOXFORHORSES@GMAIL.COM



-3-



-4-



-5-



-6-

THE INSTRUCTIONS ON THE BACK ARE YOUR RECIPE FOR SUCCESS. TAKE TIME TO READ THEM SO THAT YOU AND YOUR HORSE WILL ENJOY THE BENEFITS OF THE SOCKS.

DO NOT ride or exercise your horse while it is wearing Silver Whinnys without the use of a sport boot to help the socks stay in place. Bell boots will not work. Check them often until you are sure the sport boot is holding them in place.

Sox For Horses, Inc. disclaims any and all liability for damages, injuries or death To persons, property or animals resulting from failure to follow instructions, or use Silver Whinnys® properly, or for injuries resulting from efforts to place Silver Whinnys® on an equine.