


## Custom Silver Whinnys® for the Draft, Draft-Cross, and Elephants

We can bandage an elephant leg if given the chance. Sock knitting machines are limited to the diameter of a sock that can be knit by the cylinder size, but that doesn't stop us from creating socks to fit large legs and go over large hooves. We provide instructions at no cost that can help you use Silver Whinnys to make your own custom socks, or we can make them for you. See [How To Build Custom Silver Whinnys® Draft Size Sox](#) or [contact us](#) to order the needed sox or for our sewing services.



*Our custom draft sock.*

See: [Circulation, Not Compression, Critical to Cellulitis, Lymphangitis Relief](#)

[“Four Stages of Healing & Best Use of Silver Whinnys®”](#)

Also See [Sock Support](#)

“Hook and Loop Velcro Silver Whinnys®/article #17”

“About Compression and Silver Whinnys®/ article #06”

And If you wish to make your own hook & loop closure Silver Whinnys go to:

["How to Make Hook and Loop Silver Whinnys"](#)

